

## AMPure XP Magnetic Beads PCR Purification Process

---

We use AMPure XP Magnetic Beads to purify DNA and PCR products for genomic applications. This will remove unused dNTPs, enzymes, etc. The protocol below is modified from the product manual. The beads can be purchased from Beckman Coulter:

<https://www.beckman.com/reagents/genomic/cleanup-and-size-selection/pcr>

---

1. Make enough new 70% ethanol so you have at least 400  $\mu$ l per sample to be purified.  
*-80% ethanol can also be used as indicated in Nanopore protocol.*
  2. Shake the AMPure XP Magnetic Beads very well before use and shake again before adding them to the DNA. The color of the mixture should appear homogeneous before application.  
*-The beads fall out of suspension very quickly so shaking them well immediately before application is critical.*
  3. Add 1.8  $\mu$ l of AMPure XP beads per 1  $\mu$ l of sample and let the mixture rest for 5 minutes at room temperature.  
*-For 30  $\mu$ l of DNA, add 54  $\mu$ l of beads (30 X 1.8).*
  4. Put the mixture on a magnetic rack and let it sit for at least 2 minutes until the beads are concentrated against the magnet. The liquid should appear clear.
  5. Remove the cleared solution from the tube and discard. Leave 5  $\mu$ l of liquid in the tube so the beads are not drawn out with the supernatant. Do not disturb the beads!
  6. Dispense 200  $\mu$ l of newly made 70% ethanol alcohol to each sample and incubate at room temperature for 30 seconds. Remove the ethanol and discard.
  7. Repeat for a total of two washes.
  8. After the second wash, remove as much of the ethanol as possible without disturbing the pellet and let dry for at least 30 seconds. Do not let the pellet dry to the point of cracking.
  9. Add 30  $\mu$ l of molecular biology grade water (or needed volume) and incubate for two minutes.
  10. Pellet the beads on the magnetic rack until the liquid is clear and colorless. This should be at least 1 minute.
  11. Remove and retain the 30  $\mu$ l of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.
-